

Crusade

New possibility in the era of DES for bifurcation lesion.


1 *Optimized configuration and material achieved high shaft maneuverability.*

- Thin and flexible tip : **Distal shaft**
- Pliable and sturdy : **Proximal shaft**

2 *Configuration of double-layered lumens allows more flexible GW movement*

Product Specifications

Catalogue number	Maximum diameter for the compatible guidewire	Effective length	
KMF0114A	0.014inch (0.36mm)	1400mm	



 **KANEKA
MEDICAL
PRODUCTS**

CE 0123

KANEKA CORPORATION : 3-2-4, Nakanoshima, Kita-ku, Osaka, Japan : Phone +81-6-6226-5285 / Fax +81-6- 6226-5143