




ASAHI Meister 16

ASAHI peripheral vascular guide wire


Structure & ordering information



Product	Catalog No.	Usable length	Shaft O.D.	Radiopaque length	Polymer jacket length	Hydrophilic coating length	Tip shape
ASAHI Meister16	WAMS-165-1645	165 _{cm}	0.41 _{mm} (0.016 _{inch})	5 _{cm}	160 _{cm}	160 _{cm}	Angled 45°
	WAMS-180-1645	180 _{cm}	0.41 _{mm} (0.016 _{inch})	5 _{cm}	175 _{cm}	175 _{cm}	Angled 45°
	WAMS-165-16WA	165 _{cm}	0.41 _{mm} (0.016 _{inch})	5 _{cm}	160 _{cm}	160 _{cm}	Angled double
	WAMS-180-16WA	180 _{cm}	0.41 _{mm} (0.016 _{inch})	5 _{cm}	175 _{cm}	175 _{cm}	Angled double
	WAMS-165-16ST	165 _{cm}	0.41 _{mm} (0.016 _{inch})	5 _{cm}	160 _{cm}	160 _{cm}	Round curve
	WAMS-180-16ST	180 _{cm}	0.41 _{mm} (0.016 _{inch})	5 _{cm}	175 _{cm}	175 _{cm}	Round curve



Angled 45°



Angled double



Round curve



Not all shapes and sizes are available in all markets. Please check availability with your sales representative.
Coated with SLIP-COAT® coating.

© 2017 ASAHI INTECC CO., LTD.
"ASAHI" and "ASAHI Meister" are trademarks or registered trademarks of ASAHI INTECC CO., LTD. in Japan and other countries.

Manufactured by:

Distributed by:



ASAHI INTECC CO.,LTD.

3-100 Akatsuki-cho, Seto, Aichi 489-0071 Japan

Contact phone number: +81-561-86-9101

<http://www.asahi-intecc.com>