
Y-Connector


OKAY II™

Y-Connector Set



Connecting
design,
function and
performance

Features



Hemostatic valve

The hemostatic valve can be opened by pushing the opener, and closed by releasing it.

Once the hemostatic valve has been opened, a clockwise quarter turn locks the valve in the open position.

Fixed valve

Valve to fix devices and prevent loss of pressure.

Large lumen

Nominal ID: 10F, 3.33mm

Allows for use in combination with various products for a variety of PCI strategies.

Superior design

Facilitates single-handed operation.

Ease-of-use

Hemostatic valve features a locking mechanism.

Minimizes blood loss

Two valves (hemostatic and fixed) greatly reduce blood loss during the procedure.

Accessories

Y-Connector



Torque Device



Inserter



A set (YOK0A)

Y-connector, Inserter, Torque Device

■ Distributed by

ASAHI INTECC

ASAHI INTECC CO., LTD., EUROPE OFFICE
Strawinskylaan 527 (WTC Tower B - 5th floor) 1077 XX
Amsterdam, The Netherlands
TEL: 31-(0)20-794-0642 FAX: 31-(0)20-794-0641

GOODTEC Y-CONNECTOR SET

Manufacturer: Goodman

▶ Please consult your sales representative regarding stock availability.

▶ OKAY II is a trademark of Goodman Co., Ltd.

▶ Please refer to the Instructions For Use for further details regarding product use and specifications.

CE
2797

OKAY II-1906-A1