

Technical details

Introducer Sheath

Material – outer	Stainless steel reinforced polymer		
Material – inner	PTFE liner		
Flexibility	Regular & Firm		
Sheath coating	Hydrophobic		
Sheath tip	Smooth taper with radiopaque platinum iridium marker band		
French size	4, 5, 6, 7, 8 F		
Length	15, 45, 60, 90 cm		
Sidearm/hemostatic valve	Detachable, with color coded 3-way-stopcock		

Ordering information

Diameter (French)	Length (cm)	Catalogue Nur Regular	mber Firm new	Inner Diameter Minimum (mm/inch)	Dilator Guidewire Compatibility max.	Units/ Case
4 F	15 cm	FG-04015		1.60 mm / 0.063"	0.035" & 0.018"	5
4 F	45 cm	FG-04045		1.60 mm / 0.063"	0.035"	5
4 F	60 cm	FG-04060		1.60 mm / 0.063"	0.035"	5
4 F	90 cm	FG-04090		1.60 mm / 0.063"	0.035"	5
5 F	15 cm	FG-05015		1.92 mm / 0.075"	0.035" & 0.018"	5
5 F	45 cm	FG-05045		1.92 mm / 0.075"	0.035"	5
5 F	60 cm	FG-05060		1.92 mm / 0.075"	0.035"	5
5 F	90 cm	FG-05090		1.92 mm / 0.075"	0.035"	5
6 F	15 cm	FG-06015	FG-06015-F	2.25 mm / 0.088"	0.035" & 0.018"	5
6 F	45 cm	FG-06045	FG-06045-F	2.25 mm / 0.088"	0.035"	5
6 F	60 cm	FG-06060	FG-06060-F	2.25 mm / 0.088"	0.035"	5
6 F	90 cm	FG-06090	FG-06090-F	2.25 mm / 0.088"	0.035"	5
7 F	45 cm	FG-07045		2.55 mm / 0.100"	0.038"	5
7 F	60 cm	FG-07060		2.55 mm / 0.100"	0.038"	5
7 F	90 cm	FG-07090		2.55 mm / 0.100"	0.038"	5
8 F	45 cm	FG-08045		2.88 mm / 0.113"	0.038"	5
8 F	60 cm	FG-08060		2.88 mm / 0.113"	0.038"	5
8 F	90 cm	FG-08090		2.88 mm / 0.113"	0.038"	5

15 cm versions contain two dilators. Additional sizes to become available.








Kardia Srl Via Cormons, 18 20151 Milano (MI) www.kardia.it

DuraSheath

Kink-resistant introducer sheath

Designed for exceptional kink-resistant performance and non-ovaling effect.



Supporting complex vascular access

DuraSheath

Kink-resistant introducer sheath

From easy insertion to smooth navigation to superb visibility: DuraSheath is designed to combine optimal performance with excellent pushability and trackability.

While maintaining the inner diameter, DuraSheath allows for smooth device passage and has an exceptionally small kink-radius and non-ovaling effect. The detachable hemostatic valve provides flexibility and versatility, giving you the ability to focus on the intervention.



New optional shaft material optimal backup support transitioning to softer distal zone for safe navigation.



Dilator with snap-fit secures the dilator for safe insertion.



Detachable hemostatic valve provides you options.



properties, no matter if wet or dry.

Rethink the basics

Passionate about bringing innovative ideas and medical device solutions to life, we are a team building on nearly 20 years of experience in the design and manufacturing of medical device solutions for partners around the world.

Even small changes in product design and quality can make a clear difference – making procedures easier and simpler with better outcomes.

Continually striving for solutions to the most basic problems, we don't go back to basics, we rethink the basics.



Outer diameter

with smooth surface for low friction navigation.



Inner diameter

smooth device delivery and movement.



Hydrophobic coating

allowing for effortless & durable gliding

designed for smooth and easy insertion.

Platinum iridium marker band for excellent radiopacity.